


## List of ELISA Kits and Reagents for Food Safety and Drug Residues

ADI has highly specialized Drug residue testing ELISA kits for animal, human, or aquatic samples.

(See Details at the website) [http://4adi.com/commerce/catalog/spcategory.jsp?category\\_id=2739](http://4adi.com/commerce/catalog/spcategory.jsp?category_id=2739)

Product Name	Cat#	# tests Assay/ Range (ppb)	Samples	Detection Limit (ppb)	Sample Processing Time	Assay time
<b>Clenbuterol ELISA Test Kit</b>	DE-100010	96 tests 0.1-8.1 ppb	Pork/Liver	0.1	50 min	80 min
			Serum	0.1	10 min	80 min
			Urine	0.1	10 min	80 min
			Feed	50	30 min	80 min
			Meat	0.2	20 min	80 min
<b>Ractopamine ELISA Test Kit</b>	DE-100020	96 tests 0.1-8.1 ppb	Liver	0.4	20 min	80 min
			Urine	0.2	10 min	80 min
			Feed	15	30 min	50 min
<b>Salbutamal ELISA Test Kit</b>	DE-100030	96 tests 0.5-40.5 ppb	Urine	0.5	10 min	80 min
			Tissue	0.5	20 min	80 min
			Feed	50	30 min	80 min
<b>Chloramphenicol ELISA Test Kit</b>	DE-100040	96 tests 0.05-4.05 ppb	Animal tissue	0.05	20 min	80 min
			Aquatic	0.05	20 min	80 min
			Honey	0.15	20 min	80 min
			Intestine	0.1	20 min	80 min
			Urine	0.1	150 min	80 min
			Egg	0.1	30 min	80 min
			Milk	0.05	30 min	80 min
			Serum	0.1	15 min	80 min
<b>Florfenicol ELISA Test Kit</b>	DE-100050	96 tests 0.5-40.5 ppb	Animal tissue	0.5	30 min	80 min
			Aquatic	0.5	30 min	80 min
			Honey	1.5	30 min	80 min
			Fish, Shrimp	0.3	120 min	80 min
<b>Nitrofurans (AMOZ) ELISA Test Kit</b>	DE-100060	96 tests 0.1-8.1 ppb	Honey	0.2	120 min	80 min
			Chicken/Liver	0.2	120 min	80 min
<b>Nitrofurans (AHD) ELISA Test Kit</b>	DE-100070	96 tests 0.05-4.05 ppb	Fish, Shrimp	0.3	16 h	80 min
			Honey	0.2	16 h	80 min
			Chicken/Liver	0.2	16 h	80 min
<b>Nitrofurans (AOZ) ELISA Test Kit</b>	DE-100080	96 tests 0.05-4.05 ppb	Fish, Shrimp	0.2	16 h	80 min
			Honey	0.1	16 h	80 min
			Chicken/Liver	0.1	16 h	80 min
<b>Sulfonamides residues (SAs) ELISA Test Kit</b>	DE-100090	96 tests 1-81 ppb	Chicken/Liver	1	20 min	90 min
			Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
			Milk	20	10 min	90 min
<b>Sulfamethazine (SM2) ELISA Test Kit</b>	DE-100100	96 tests 1-81 ppb	Chicken/Liver	1	20 min	90 min
			Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
			Milk	20	10 min	90 min
<b>Sulfameth- oxydiazine (SMD) ELISA Test Kit</b>	DE-100110	96 tests 1-81 ppb	Chicken/Liver	1	20 min	90 min
			Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
<b>Sulfaquin- oxaline ELISA Test Kit</b>	DE-100240	96 tests 1-81 ppb	Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
			Milk	20	10 min	90 min



## List of ELISA Kits and Reagents for Food Safety and Drug Residues

ADI has highly specialized Drug residue testing ELISA kits for animal, human, or aquatic samples.

(See Details at the website) [http://4adi.com/commerce/catalog/spcategory.jsp?category\\_id=2739](http://4adi.com/commerce/catalog/spcategory.jsp?category_id=2739)

Product Name	Cat#	# tests Assay/ Range (ppb)	Samples	Detection Limit (ppb)	Sample Processing Time	Assay time
<b>Quinolones (QNS) ELISA Test Kit</b>	DE-100120	96 tests 1-81 ppb	Pork/Liver	3	50 min	90 min
			Chicken/Liver	3	50 min	90 min
			Shrimp, Fish	3	20 min	90 min
			Serum	3	10 min	90 min
			Honey	6	10 min	80 min
<b>Enrofloxacin ELISA Test Kit</b>	DE-100130	96 tests 0.25-20.25 ppb	Pork/Liver	1	50 min	80 min
			Chicken/ liver	1	50 min	80 min
			Shrimp, Fish	1	20 min	80 min
			Serum	1	10 min	80 min
			Honey	2	10 min	80 min
<b>Ampicillin ELISA Test Kit</b>	DE-100140	96 tests 0.1-8.1 ppb	Pork/Liver	2	15 min	80 min
			Chicken (meat)	2	15 min	80 min
			Duck (meat)	2	15 min	80 min
			Shrimp, Fish	2	15 min	80 min
			honey, milk	2	10 min	80 min
<b>Benzyl penicillin ELISA Test Kit</b>	DE-100150	96 tests 0.1-8.1 ppb	Pork/liver	2	15 min	80 min
			Chicken(meat, liver)	2	15 min	80 min
			Duck(meat, liver)	2	15 min	80 min
			Shrimp, Fish	2	15 min	80 min
			Honey, Milk	2	10 min	80 min
<b>Tylosin ELISA Test Kit</b>	DE-100160	96 tests 1.5-121.5 ppb	Meat	3	20 min	60 min
			Liver	3	20 min	60 min
			Honey/Egg	1.5	10 min	60 min
<b>Trenbolone ELISA Test Kit</b>	DE-100170	96 tests 0.1-8.1 ppb	Animal tissue	5	30 min	80 min
			Aquatic	5	30 min	80 min
			Urine	0.5	30 min	80 min
<b>Diazepam ELISA Test Kit</b>	DE-100180-	96 tests 0.1-8.1 ppb	Tissue	1	30 min	80 min
			Urine	1	30 min	80 min
			Feed	20	30 min	80 min
<b>Diethylstilbestrol (DES) ELISA Test Kit</b>	DE-100190	96 tests 0.1-8.1 ppb	Fish, Shrimp	0.2	20 min	80 min
			Liver, Meat	2	20 min	80 min
			Feed	20	20 min	80 min
			Urine	0.6	20 min	80 min
<b>Gentamicin ELISA Test Kit</b>	DE-100200	96 tests 0.1-8.1 ppb	Chicken/Liver	4	40 min	80 min
<b>Streptomycin ELISA Test Kit</b>	DE-100210	96 tests 0.1-25.6 ppb	Chicken/Liver	20	40 min	60 min
			Honey	10	40 min	60 min
			Milk	20	5 min	60 min
<b>Nitrofurantoin (SEM) ELISA Test Kit</b>	DE-100070	96 tests 0.1-8.1 ppb	Honey	0.2	16 h	80 min
			Fish, Shrimp	0.3	16 h	80 min
			Chicken/Liver	0.2	16 h	80 min
			Fish, Shrimp	0.3	16 h	80 min
			Chicken/Liver	0.2	16 h	80 min
<b>Olaquinox ELISA Test Kit</b>	DE-100230	96 tests 1-81 ppb	Tissue	1	40 min	80 min
			Feed	100	45 min	80 min

Rev. 130118A

